

Now available with reflective backing!
Try our new Reflective DryTips®

Video and instructions for use available
at www.microbrush.com/drytips

Ultimate Saliva Control

Far outlasts cotton rolls and other absorbent pads, and doesn't lose moisture when handled.

The DryTips® advantage
over traditional cotton rolls:

- Complete coverage of parotid orifice and entire buccal mucosa with polyethylene film for cutting protection
- Facilitates work far back in the oral cavity
- Flexible; accommodates to cheek movements
- Comfortably stays in place, yet is easily removed with water spray
- Reflective coating for increased visibility and to aid in retraction



NEW! Reflective coating
now available.

Ideal for:

- Impressions
- Occlusal Adjusting
- Orthodontics
- Paint-on techniques
- Composites
- Cementation
- Fissure sealing
- Fluoride applications

DryTips® Ordering Information:

ITEM NUMBER	SIZE	QUANTITY	COLOR
291542	Large	50	White
291543	Small	50	White
291726	Large	50	Silver
291727	Small	50	Silver

Order today through
your preferred dealer.

For questions, call
Microbrush® International
at (866) 866-8698.

Your Choice of 3 Great Offers!

BUY 2

Two 4-packs of the
Microbrush® Plus
Series or Tube Series
Applicators



GET 1

One Box of New!
Reflective DryTips
for no charge



BUY 1

One Refill Pack
of Coated
or Uncoated
ContactPro® SmartBand
Matrix Bands



GET 1

One Refill Pack of the
same ContactPro
SmartBand Matrix Bands
for no charge*



*See side 2 for product details

ContactPro SmartBands can
be used with any matrix ring!

MUST WRITE Promo Code:
Q417 on Invoice

BUY 3

Three Boxes
of DryTips® or
Reflective DryTips®



GET 1

One Box of
Reflective DryTips
for no charge



Full Name		Date of Purchase	
Street Address			
City	State	Zip	
E-Mail	Telephone	Fax	

*To redeem this offer, send a copy of your qualifying invoice to: freegoods@microbrush.com or fax it to 888-771-7573. Purchases must be shown on one invoice. Offer limited to one per office and must be received no later than one month after the promotion expires. Offer good 10/1/2017-12/31/2017.

CONTACT PRO[®]

MATRIX SYSTEM

See the ContactPro[®] matrix system in action at microbrush.com/contactpro



Simple Technique for Outstanding Results

The reliable, consistent and easy to use system that provides solutions for Class II Restorations.

1 Place the Matrix Band



Select a band that matches the height of the tooth and place it with the concave edge toward the occlusal margin using the SmartTab at the top of the band.

2 Place the Wedge



Place a finger on the band to prevent dislodgment and firmly insert wedge to seal the gingival margin.

3 Place the Ring




Open the ring with your ContactPro[®] Ring Forceps and place it directly over the wedge, seating it as far gingivally as possible.

4 Voila!



Simple technique for outstanding results.

CONTACT PRO SmartBand



ContactPro[®] SmartBand Matrix Bands
Refill Pack Ordering Information:

ITEM NUMBER	SIZE	QUANTITY	COLOR
291704	3.5mm	100ct	White
291591	4.5mm	100ct	Green
291593	5.5mm	100ct	Blue
291705	6.5mm	25ct	Purple
291706	7.5mm	25ct	Gray
291707	3.5mm	100ct	Silver
291590	4.5mm	100ct	Silver
291592	5.5mm	100ct	Silver
291708	6.5mm	25ct	Silver
291709	7.5mm	25ct	Silver

Order today through your preferred dealer.
For questions, call Microbrush[®] International at (866) 866-8698.